

# Ocean-analyze Guide

Author: Kim Miikki

Date: 2.11.2021

## 1 Introduction

The spectral data parsed with *ocean2csv.py* can be analyzed with the *ocean-analyze.py*. This program will create ordinary or transmission spectra as plots and text data.

## 2 System Requirements

Operating System: ALL

Python 3 with NumPy, Matplotlib, Pandas and scikit-learn.

## 3 Spectral Methods

This program supports two different spectral methods: normal and transmission. The first method creates mean spectra plots, where y axis is Intensity in counts. The second method creates spectral mean transmission plots. The mean spectrum is calculated from accepted scans in its scan range (= scans/sample). The acceptance is controlled with an area under curve threshold. The threshold is at default (2/3 of highest area or 2/3 of highest area mean of 1...N-1 areas). When the area falls under the mean area, the spectral data is discarded within other smaller-area data in the same scan series. An area plot is created when more than one spectrum is analyzed.

## 4 Program Usage

The program arguments are listed in Table 1.

Argument	Description
-h, --help	show this help message and exit
-s S	scans/spectrum
-t T	threshold value: 0.01-0.99
-d	display images on screen
-i	ignore threshold
-g	enable grid
-y	auto y axis
-ts	transmission spectra

Table 1. Program arguments.

At least scans/spectrum is asked if running the program without any arguments. The first spectrum is analyzed, and if trying to analyze transmission spectra as ordinary spectra, a data type question is asked if the transmission degree is at least 5 %. The grid can be enabled with the argument -g and auto y axis with the argument -y. When analyzing transmission spectra, the argument -ts should be used.

## 5 Use Cases

The usage of this program is demonstrated in the following subsections.

### 5.1 Analyzing Normal Spectra Files

Spectra parsed with *ocean2csv.py* were analyzed with *ocean-analyze.py*:

```
$ ocean-analyze.py
```

```
Ocean spectrometer CSV spectra analyzer
```

```
File format: CSV with a header line
```

```
Current directory:
```

```
/home/pi/python/20211102-ocean/data
```

```
Select scans/spectrum: (1...35, Default=1: <Enter>): 5
```

```
Processing:
```

```
00-20211021_FLMT044961__0__16.csv
```

```
01-20211021_FLMT044961__1__17.csv
```

```
02-20211021_FLMT044961__2__18.csv
```

```
03-20211021_FLMT044961__3__19.csv
```

```
04-20211021_FLMT044961__4__20.csv
```

```
00-spectrum-mean_(5#5).csv
```

```
...
```

```
30-20211021_FLMT044961__30__46.csv
```

```
31-20211021_FLMT044961__31__47.csv
```

```
32-20211021_FLMT044961__32__48.csv
```

```
33-20211021_FLMT044961__33__49.csv
```

```
34-20211021_FLMT044961__34__50.csv
```

```
30-spectrum-mean_(5#5).csv
```

```
Saving area files
```

The scans value was 5, therefore a set of five spectra is analyzed and compared with the area under curve. An analysis subdirectory is created for analysis files. In this case following files were created:

```
$ ls
```

```
'00-spectrum-mean_(5#5).csv' '20-spectrum-mean_(5#5).csv'
```

```
'00-spectrum-mean_(5#5).png' '20-spectrum-mean_(5#5).png'
```

```
'05-spectrum-mean_(5#5).csv' '25-spectrum-mean_(4#5).csv'
```

```
'05-spectrum-mean_(5#5).png' '25-spectrum-mean_(4#5).png'
```

```
'10-spectrum-mean_(5#5).csv' '30-spectrum-mean_(5#5).csv'
```

```
'10-spectrum-mean_(5#5).png' '30-spectrum-mean_(5#5).png'
```

```
'15-spectrum-mean_(5#5).csv' spectra-areas.csv
```

```
'15-spectrum-mean_(5#5).png' spectra-areas.png
```

The mean spectrum per scan series is stepped by 5. The first number correspond to the scan set. The PNG files contains the plots and CSV the data. The first spectrum is shown in Fig. 1.




Figure 1. Blue LED on the left and its spectrum on the right.

The spectrum data is saved in a csv file:

```
$ head 00-spectrum-mean_(5#5\).csv
Wavelength (nm),Intensity (counts)
179.018,-271.09
179.241,-271.09
179.464,-271.09
179.687,-249.606
179.91,-12.238
180.133,-28.44
180.356,18.048
180.578,-28.438
180.801,27.912
```

The 5#5 indicate that 5 spectra of 5 are accepted to the mean spectrum. The area curve for all scans are shown in Figure 2.




Figure 2. Spectra areas when the number of scans/series is 5.

## 5.2 Analyzing an IR-pass Filter in Transmission Mode

The transmittance of an infrared pass filter was measured with an Ocean spectrometer, the data was parsed with *ocean2csv.py* and finally the data was analyzed and plotted with *ocean-analyze.py*:

```
$ ocean-analyze.py -g -ts -s 5
```

Ocean spectrometer CSV spectra analyzer

File format: CSV with a header line

Current directory:

/home/pi/python/20211102-ocean/bw/data

Processing:

0-20211021\_Transmission\_\_35\_\_35.csv

1-20211021\_Transmission\_\_36\_\_36.csv

2-20211021\_Transmission\_\_37\_\_37.csv

3-20211021\_Transmission\_\_38\_\_38.csv

4-20211021\_Transmission\_\_39\_\_39.csv

0-transmission-spectrum-mean\_(5#5).csv

qt5ct: using qt5ct plugin

Saving area files

The mean transmission spectrum is shown in Figure 3.




Figure 3. B+W 092 (DARK RED 20-40X) ø 55 mm filter transmission spectrum.